

Protein Extraction



Medicago truncatula



Shinorhizobium medicae
root nodules

Plant Harvest:

drought stress



Control



Liquid-N2!

Grinding and solubilizing:

Urea, Thiourea, CHAPS, Triton X-100, DTT, Ampholytes, Trizma-base, Trizma-HCL



Grinding

Centrifugation



Protein extract

various extraction buffers and detergents exists!

For private use only!

Dept. Molecular Systems Biology

PP Proteomics (Wienkoop)

Protein Extraction, Desalting, Digestion



Nur zum persönlichen Gebrauch!

Protein Extraction



Medicago truncatula



Sinorhizobium medicae
root nodules

Plant Harvest:

drought stress



Control



Medicago truncatula (*Rhizobium* inoculated) : **leaves**

Group 1) wt

Group 2) *dnf5-2* (N-fixation deficient)

Group 3) NF11301 (ferritin F3)

Group 4) NF9644 (ferritin F2)




each: **control** and **drought stress** (7 days)

GLP



Skin flakes

Dandruff



GLP

clean your bench!



use the gloves!

For private use only!

Dept. Molecular Systems Biology

PP Proteomics (Wienkoop)

Protein Extraction, Desalting, Digestion



Nur zum persönlichen Gebrauch!

GLP

Handle with care!!!



NEVER !!!



For private use only!

Dept. Molecular Systems Biology

PP Proteomics (Wienkoop)

Protein Extraction, Desalting, Digestion



Nur zum persönlichen Gebrauch!

Safety Instructions!!!



Lena Fragner



For private use only!

Dept. Molecular Systems Biology

PP Proteomics (Wienkoop)

Protein Extraction, Desalting, Digestion



Nur zum persönlichen Gebrauch!

Different Extraction strategies

- ❖ **Volatile buffers preferable (e.g. ammonium bicarbonate)**
 - can be removed by vacuum centrifugation or desalting
 - reconstituted in aqueous organic solvent + acid
- ❖ **Soluble samples**
 - 50 % aqueous MeOH or CH₃CN + 1-5 % formic/acetic acid
- ❖ **Insoluble samples (membrane proteins)**
 - chloroform, higher percentage of acid

There is no „one and only“ extraction!

Protein Extraction

- ★ Urea for protein solubilization; does not bind to proteins – can cause modification at temperatures above 37°C!
- ★ PMSF (Phenylmethylsulfonylfluorid) = Protease inhibitor
- ★ Acetone Precipitation for desalting and removal of hydrophobic and other MS-unsuitable buffer components
- ★ β -Mercaptoethanol like DTT for irreversible disulfide reduction

Protein Extraction

**Gel-free
vs.
Gel-based**

Shotgun-MS Analysis

For private use only!

Dept. Molecular Systems Biology

PP Proteomics (Wienkoop)

Protein Extraction, Desalting, Digestion




Nur zum persönlichen Gebrauch!


Protein Digestion

For **gel-free** Shotgun-MS Analysis

In Solution
Protein digestion:




Protein



Peptides


Peptide desalting:



Peptides in Urea
(Salt) Buffer

SPEC18

Solid Phase Extraction
using
Reversed Phase



For private use only!

Dept. Molecular Systems Biology

PP Proteomics (Wienkoop)

Protein Extraction, Desalting, Digestion



Nur zum persönlichen Gebrauch!

Protein Extraction

- ❖ **Removal of non-volatile buffers i.e. NaCl, phosphate, etc**
 - reverse phase HPLC (off-line)
 - desalting cartridges
 - centrifugal filtration
 - passive dialysis
 - interface chromatographic or desalting cartridges with ESI source (LC-MS)

- ❖ **Removal of detergents or surfactant**
 - they quench and obscure ESI response (e.g. SDS from gel)
 - difficult to remove from protein samples, SCX
 - best to avoid if sample is for ESI analysis

Peptide Desalting


Sample purification

- a) activate – MeOH
- b) equilibrate – solution A
- c) load you sample
- d) wash – solution A
- b) elute – solution B


solution A:
0.5% acetic acid

solution B:
0.5% acetic acid
80% acetonitrile

How it works ?



Sample preparation workflow



Lysate samples are prepared from biological specimens or cultured cells by a customized protocol that may include cell lysis, subcellular fractionation, depletion of high abundant proteins, enrichment of target proteins, dialysis and desalting. In-solution digestion entails irreversibly breaking disulfide bonds via reduction and alkylation followed by protein digestion into peptide fragments. An alternative approach is to first resolve proteins by 1D or 2D electrophoresis (1DE or 2DE, respectively) and then collect gel slices that contain the desired band(s). The proteins in these gel plugs are then reduced, alkylated and digested in situ. After peptides are extracted from the gel matrix, the peptides are enriched and salts and detergents removed and prepared for MS analysis. While this diagram encompasses the most common steps of sample preparation, the protocol must be tailored to specific samples and MS techniques for optimal results.

For private use only!

Dept. Molecular Systems Biology

PP Proteomics (Wienkoop)

Protein Extraction, Desalting, Digestion




Nur zum persönlichen Gebrauch!

Protein Extraction


For **gel-based** Shotgun-MS Analysis

In Gel
Protein digestion:

SDS- PAGE



2D-(IF-SDS) Gel Electrophoresis



or

- 1) Cutting gel band/spot
- 2) Wash gel pieces
- 3) In-Gel Tryptic Digestion
- 4) Peptide extraction from gel
- 5) Desalting with STAGE Tips
- 6) MS-Analysis

Sample preparation workflow

digest

Cut the band out



Wash

Volatile buffers - ammonium bicarbonate!

Reduction

DTT

Alkylation

Iodoacetamide

Digestion

Trypsin or LysC

Extraction

Acetonitrile

Purification

STAGE Tip



MS

For private use only!

Dept. Molecular Systems Biology

PP Proteomics (Wienkoop)

Protein Extraction, Desalting, Digestion





Nur zum persönlichen Gebrauch!

Reduction and Alkylation

Digest


Reduction



DTT to reduce SH-bridges

Alkylation

Carbamidomethylation of Cysteine



Peptide Desalting


SPE C18 (Solid Phase Extraction)



STop And Go Extraction Tip (StageTip) vs. 96-well-SPEC

STAGE tip is a reversed-phase column providing:


- a) Desalting
- b) Filtration
- c) Concentration



From www.mmm.co.jp

Sample preparation workflow

LC/MS/MS Shotgun Analysis



➔ Screening for marker enzymes

Good Laboratory Practice

**BEFORE YOU START !!!
MAKE SURE !!!
THAT YOU AVOID !!!**

For private use only!

Dept. Molecular Systems Biology

PP Proteomics (Wienkoop)

Protein Extraction, Desalting, Digestion



Nur zum persönlichen Gebrauch!